



Multistep HTC-D2B PROTAC Synthesis Enabled by Late-Stage sp^2 - sp^3 Cross-Coupling Chemistry



Rebecca Stevens^{1,2}, Afjal Miah¹, Glenn Burley²

1) GSK, Medicines Research Centre, Stevenage, UK

2) Department of Pure and Applied Chemistry, University of Strathclyde, Glasgow, UK

Develop late-stage $C(sp^2)$ - $C(sp^3)$ cross coupling conditions



Apply conditions to generate PROTAC intermediates



Synthesise diverse PROTAC library using D2B platform




Identify novel PROTAC hits

AIMS


- Identify $C(sp^2)$ - $C(sp^3)$ cross-coupling conditions suitable for use in a multistep High-Throughput-Chemistry Direct-to-Biology (HTC-D2B) workflow to expand the GSK PROTAC platform capabilities.
- Utilise new conditions to synthesise diverse libraries of PROTACs, exploring wider chemical space and physchem properties.

LATE-STAGE FUNCTIONALISATION





Conversion
High
Low

- High-throughput experimentation aided identification of suitable conditions.
- Screened catalysts and % loading, bases, solvents, temperature, reaction time.




- Substrate scope investigated in plate format - conditions suitable for 3 different fully-functionalised POI binders and a wide range of linkers.
- Reactions scaled up with BRD4 binder in moderate to good yields:




Results inform future PROTAC design

Structure-activity relationship (SAR) analysis




Computational data analysis

Late-stage functionalisation step to synthesise structurally diverse intermediates



Automated PROTAC synthesis and LCMS analysis in 1536- and 384-well plates

High-Throughput Chemistry Direct-to-Biology
HTC-D2B PLATFORM





Degradation data for thousands of novel PROTACs rapidly obtained

- 700 reactions carried out with chemistry success rate of ~70%, to synthesise 500 PROTACs in one month - assay to be run in the next month. Platform fully validated for a first set of 74 PROTACs, biological data below.


IDENTIFICATION OF NOVEL BRD4 HITS

- First set of 74 PROTACs identified many highly potent compounds, several with desirable physchem properties for oral bioavailability.
- Identified novel hits for BRD4, employing both cereblon and VHL.
- 4 hits resynthesised and purified as controls - excellent correlation between D2B and purified samples.

DATA FOR D2B VS. PURIFIED SAMPLES



*pDC₅₀ of 9.77 = limit of assay (rounded to 1 dp in bar chart); colours represent two series of intermediates; n=1 for all assay data presented.



CONCLUSION

- High-throughput experimentation enabled chemistry optimisation of a late-stage $C(sp^2)$ - $C(sp^3)$ coupling.
- 500 PROTACs have been synthesised in one month, with 74 novel BRD4 PROTACs already tested in HiBit assay.
- D2B identified many highly potent hits, with purified compounds showing close correlation and no cytotoxicity.

ACKNOWLEDGEMENTS

Thanks to Diana Klimaszewska, Justyna Macina and Peter Stacey from Protein Degradation for running the first set of biological assays, and for their training in the relevant biology for future assays.